

FBXO6 antibody



Catalog Number: 110557

Product name

FBXO6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FBXO6 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, Transfected HEK-293 cells. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of FBXO6 recombinant protein (Accession Number: NM_018438). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110557(FBXO6 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110557(FBXO6 antibody) at dilution of 1:50 (under 10x lens)