

Anti-CLTCL1 antibody



Catalog Number: 176557

Product name

Anti-CLTCL1 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit monoclonal antibody to CLTCL1

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*TBS (pH7.4), 0.05% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2,000

IHC-P: 1:50-1:200

Validations




Fig1: Western blot analysis of Clathrin heavy chain on Hela (1) and PC-12 (2) cell lysate using anti-Clathrin heavy chain antibody at 1/1,000 dilution.