

# UGP2 antibody

Catalog Number: 116557

## Product name

UGP2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

UGP2 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive WB detected in HepG2 cells, HEK-293 cells, HeLa cells, mouse liver tissue. Observed molecular weight by Western-blot: 50-56kd

## Preparation

This antibody was obtained by immunization of UGP2 recombinant protein (Accession Number: NM\_001377525). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116557(UGP2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:116557(UGP2 antibody) at dilution of 1:100 (under 25x lens)