

N4-Methylcytidine Polyclonal Antibody



Catalog Number: 162957

Product name

Rabbit IgG

N4-Methylcytidine Polyclonal Antibody

Applications**Specificity**

WB, IF

Human

Dilutions**Antibody description**

WB 1:500 - 1:2000

Polyclonal antibody to N4-Methylcytidine


IF 1:50 - 1:200

Preparation

Antigen: A synthetic peptide of human 4mC

Validations**Formulation**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Immunofluorescence - N4-Methylcytidine
Polyclonal Antibody

Clonality

Polyclonal

Immunofluorescence analysis of HeLa cells using
4mC antibody at dilution of 1:100 (40x lens). Blue:
DAPI for nuclear staining.

Ig Type