

# ZNF195 antibody



Catalog Number: 116957

## Product name

ZNF195 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ZNF195 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue, human heart tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 64 kDa

## Preparation

This antibody was obtained by immunization of ZNF195 recombinant protein (Accession Number: NM\_001130519). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116957(ZNF195 Antibody) at dilution of 1:1000