

# DNA-PKcs antibody



Catalog Number: 109957

## Product name

DNA-PKcs antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

DNA-PKcs Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, MCF7 cells. Observed molecular weight by Western-blot: 350-460 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_006904). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109957(PRKDC antibody) at dilution of 1:300