

SEPT5 Polyclonal Antibody



Catalog Number: 164257

Product name

SEPT5 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to SEPT5

Preparation

Antigen: A synthetic Peptide of human SEPT5

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC


Dilutions

WB 1:500 - 1:1000

IHC 1:20 - 1:50

Validations

Human brain



Western blot - SEPT5 Polyclonal Antibody

Western blot analysis of extracts of human brain, using SEPT5 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.