

# SKIL antibody

Catalog Number: 115257

## Product name

SKIL antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SKIL Rabbit Polyclonal antibody. Positive WB detected in A431 cells, HepG2 cells, human skeletal muscle tissue. Positive IP detected in HepG2 cells. Observed molecular weight by Western-blot: 77 kDa

## Preparation

This antibody was obtained by immunization of SKIL recombinant protein (Accession Number: NM\_005414). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:115257(SKIL antibody) at dilution of 1:300



IP Result of anti-SKIL (IP:Catalog No:115257, 4ug; Detection:Catalog No:115257 1:500) with HepG2 cells lysate 1600ug.