

Anti-STAT5B Rabbit antibody



Catalog Number: 167257

Product name

Anti-STAT5B Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

STAT5B Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human STAT5B

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/1000

ICC/IF: 1/100

FC: 1/20

Validations



Western blot detection of STAT5B in C6 cell lysates using STAT5B Rabbit pAb(1:1000 diluted).Predicted band size:90kDa.Observed band size:90kDa.