

HDAC1 antibody



Catalog Number: 107257

Product name

HDAC1 antibody

Immunogen

[Human HDAC1 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

HDAC1 Mouse Monoclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, HepG2 cells, Jurkat cells. Observed molecular weight by Western-blot: 65kD

Preparation

This antibody was obtained by immunization of HDAC1 recombinant protein (Accession Number: NM_004964). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:1000-1:10000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107257(HDAC1 antibody) at dilution of 1:4000