

# Znt5 antibody

Catalog Number: 117257

## Product name

Znt5 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Znt5 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 84-87 kDa

## Preparation

This antibody was obtained by immunization of Znt5 recombinant protein (Accession Number: NM\_024055). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:200-1:2000

## Validations



IP Result of anti-SLC30A5 (IP:Catalog No:117257, 4ug; Detection:Catalog No:117257 1:500) with mouse liver tissue lysate 4000ug.