

Anti-SIRT6 (6C9) Mouse antibody



Catalog Number: 168257

Product name

Anti-SIRT6 (6C9) Mouse antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

SIRT6 (6C9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human SIRT6 protein expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB, ICC/IF, IP

Dilutions


WB: 1/500

;ICC: 1/100


Validations



Immunofluorescent analysis of HeLa cells fixed by 4% paraformaldehyde and using SIRT6 mouse mAb (dilution 1:100).



Immunoprecipitation analysis of HeLa cell lysates using SIRT6 mouse mAb.



Western blot analysis of extracts from Jurkat, MCF7, HCT116, K562, COS7, C2C12, 3T3 and HeLa cell lysates using SIRT6 mouse mAb (1:500 diluted). Predicted band size: 42,36kDa. Observed band size: 42,36kDa.