

# CENPO antibody



Catalog Number: 109257

---

## Product name

CENPO antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

CENPO Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HEK-293 cells, mouse kidney tissue, mouse ovary tissue. Observed molecular weight by Western-blot: 34kd

## Preparation

This antibody was obtained by immunization of CENPO recombinant protein (Accession Number: NM\_001322101). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109257(CENPO antibody) at dilution of 1:500