

IL10RB antibody



Catalog Number: 111657

Product name

IL10RB antibody

Specificity

human; other species not tested.

Antibody description

IL10RB Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, K-562 cells. Observed molecular weight by Western-blot: 37-41kd

Preparation

This antibody was obtained by immunization of IL10RB recombinant protein (Accession Number: NM_000628). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111657(IL10RB antibody) at dilution of 1:500