

MEPCE antibody

Catalog Number: 112657

Product name

MEPCE antibody

Specificity

Human, Mouse, Rat, Zebrafish; other species not tested.

Antibody description

MEPCE Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 74 kDa

Preparation

This antibody was obtained by immunization of MEPCE recombinant protein (Accession Number: NM_019606). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:112657(MEPCE antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:112657(MEPCE antibody) at dilution of 1:50 (under 10x lens)