

C14orf106 antibody



Catalog Number: 108657

Product name

C14orf106 antibody

Immunogen

[Human C14orf106 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

C14orf106 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 130 kDa

Preparation

This antibody was obtained by immunization of C14orf106 recombinant protein (Accession Number: XM_047431545). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108657(C14orf106 Antibody) at dilution of 1:1000