

PYGB antibody

Catalog Number: 114357

Product name

PYGB antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PYGB Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 97kd

Preparation

This antibody was obtained by immunization of PYGB recombinant protein (Accession Number: NM_002862). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114357(PYGB antibody) at dilution of 1:500



IP Result of anti-PYGB (IP:Catalog No:114357, 4ug; Detection:Catalog No:114357 1:500) with mouse brain tissue lysate 4000ug.