

# TREML4 antibody



Catalog Number: 116357

## Product name

TREML4 antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse; other species not tested.

## Applications

ELISA, WB

## Antibody description

TREML4 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, U-937 cells. Observed molecular weight by Western-blot: 50-55 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of TREML4 recombinant protein (Accession Number: NM\_198153). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:116357(TREML4 Antibody) at dilution of 1:500

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality