

CNPY3 antibody

Catalog Number: 109357

Product name

CNPY3 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CNPY3 Rabbit Polyclonal antibody. Positive IP detected in mouse kidney tissue. Positive WB detected in human brain tissue, Huh7 cells, Jurkat cells. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 30-40 kDa

Preparation

This antibody was obtained by immunization of CNPY3 recombinant protein (Accession Number: NM_006586). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations




WB result of anti-CNPY3 (Catalog No:109357, 1:1000) with 40ug Huh7 cells by Dr.Dan



Lindholm.


human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109357(CNPY3 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human brain using Catalog No:109357(CNPY3 antibody) at dilution of 1:50 (under 10x lens)

CNPY3 antibody

Catalog Number: 109357



IP Result of anti-CNPY3 (IP:Catalog No:109357, 3ug; Detection:Catalog No:109357 1:500) with mouse kidney tissue lysate 6000ug.