

Human ETFB/FP585 (His tag) recombinant protein



Catalog Number: 519357

General Information

Gene Name Synonym

Electron-transfer-flavoprotein beta polypeptide, MADD, beta-ETF.

Protein Construction

A DNA sequence encoding the Human ETFB/FP585 Met1-Ile255 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM HEPES, 150mM NaCl, pH 7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human ETFB/FP585 (His tag) recombinant protein