

Rat CD302/CLEC13A (His Tag) recombinant protein



Catalog Number: 504757

General Information

Protein Construction

A DNA sequence encoding the rat CD302 (Met1-His165) was expressed with a polyhistidine tag at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

(49.5+41.9) % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 21

Molecular Mass

The recombinant rat CD302 comprises 156 amino acids and predicts a molecular mass of 18.1 kDa. The apparent molecular mass of the recombinant protein is approximately 18 and 22 kDa in SDS-

PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat CD302 / CLEC13A Protein (His Tag) SDS-PAGE