

Anti-Profilin 2 / PFN2 antibody



Catalog Number: 105757

Product name

Anti-Profilin 2 / PFN2 antibody

Immunogen

[Human Profilin 2 / PFN2 \(His Tag\) recombinant protein](#)

Specificity

Human Profilin 2 / PFN2

Antibody description

Rabbit monoclonal to Profilin 2 / PFN2

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human Profilin 2 / PFN2 (rh Profilin 2 / PFN2; AAH18049.1; Met1-Phe140).

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 5-10 µg/ml

Validations



Lanes	A
Sample (whole cell lysate)	NCI-H1299
Sample Volume (µg/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	10 µg/ml
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.
Explanation	Developed using Odyssey imaging system. Predicted band size : 15 kDa Observed band size : 15 kDa

Profilin 2 / PFN2 Antibody, Rabbit MAb, Western blot