

# Anti-Profilin 2 / PFN2 antibody



Catalog Number: 105757

## Product name

Anti-Profilin 2 / PFN2 antibody

## Immunogen

[Human Profilin 2 / PFN2 \(His Tag\) recombinant protein](#)

## Specificity

Human Profilin 2 / PFN2

## Antibody description

Rabbit monoclonal to Profilin 2 / PFN2

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human Profilin 2 / PFN2 (rh Profilin 2 / PFN2; AAH18049.1; Met1-Phe140).

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 5-10 µg/ml

## Validations



| Items                      | Lanes  |
|----------------------------|--|
| Sample (whole cell lysate) | A  |
| Sample Volume (µg/lane)    | NCI-H1299<br>30  |
| Gel                        | 13% SDS-PAGE reducing gel  |
| Recommended Concentration  | 10 µg/ml   |
| Secondary Antibody         | Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.<br>Developed using Odyssey imaging system. |
| Explanation                | Predicted band size : 15 kDa<br>Observed band size : 15 kDa  |

Profilin 2 / PFN2 Antibody, Rabbit MAb, Western blot