

Anti-AVEN Rabbit antibody



Catalog Number: 166757

Product name

Anti-AVEN Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

AVEN Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human AVEN

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000-1/5000

IHC: 1/20-1/50

ICC/IF: 1/20-1/50

FC: 1/20

Validations



Western blot detection of AVEN in Rat Brain, MCF7, 3T3, HeLa cell lysates using AVEN Rabbit pAb (1:1000 diluted). Predicted band size: 39 kDa. Observed band size: 50 kDa.