

# Anti-TMEM43 Rabbit antibody



Catalog Number: 168757

## Product name

Anti-TMEM43 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

TMEM43 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human TMEM43

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions

1:1000-1:5000

1:50

## Validations



Western blot detection of TMEM43 in K562,C6 cell lysates using TMEM43 Rabbit mAb(1:1000 diluted).Predicted band size:45kDa.Observed band size:45kDa.