

CACNA1B antibody

Catalog Number: 108757

Product name

CACNA1B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CACNA1B Rabbit Polyclonal antibody. Positive IHC detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 262 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000718). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:50-1:500


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108757(CACNA1B antibody) at dilution of



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:108757(CACNA1B Antibody) at dilution of 1:100 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:108757(CACNA1B Antibody) at dilution of

CACNA1B antibody



Catalog Number: 108757

1:100 (under 40x lens). heat mediated antigen

retrieved with Tris-EDTA buffer(pH9).