

# DPM1 antibody



Catalog Number: 110067

## Product name

DPM1 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

DPM1 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas. Positive WB detected in mouse pancreas tissue, human brain tissue. Observed molecular weight by Western-blot: 27-32 kDa

## Preparation

This antibody was obtained by immunization of DPM1 recombinant protein (Accession Number: NM\_003859). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:110067(DPM1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110067(DPM1 antibody) at dilution of 1:50 (under 10x lens)