

Anti-IL18RAP antibody



Catalog Number: 142067

Product name

Anti-IL18RAP antibody

Immunogen

[Human IL18RAP \(His Tag\) recombinant protein](#)

Specificity

Human IL18RAP

Antibody description

Rabbit Polyclonal to Human IL18RAP

Preparation

Produced in rabbits immunized with purified, recombinant Human IL18RAP. IL18RAP specific IgG was purified by Human IL18RAP affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ELISA


Dilutions

WB: 5-20 µg/ml

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Human IL18RAP.

Validations



Lanes	A	B
Items		
Sample (whole cell lysate)	RAW264.7	K562
Sample Volume (µg/lane)	30	30
Gel	13% SDS-PAGE reducing gel	
Recommended Concentration	5-20 µg/ml	
Secondary Antibody	Dylight 800 labeled Antibody to Rabbit IgG (H+), at 1:5000 dilution.	
Developed using Odyssey imaging system.		
Explanation	Predicted band size: 68 kDa Observed band size: 66 kDa Additional bands at: 21 kDa, 38 kDa, 49 kDa and 50 kDa (We are unsure as to the identity of these extra bands.)	

IL18RAP Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot