

Catalog Number: 165067

Product name

ASIP Polyclonal Antibody

IHC 1:50 - 1:100

IF 1:50 - 1:100

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to ASIP

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 23-132 of human ASIP (NP_001663.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality


Polyclonal

Ig Type

Rabbit IgG


Applications

IHC, IF

Dilutions**Validations**

Immunohistochemistry - ASIP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate cancer using ASIP antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ASIP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using ASIP antibody at dilution of 1:100 (40x lens).