

Rolinsatamab reference antibody



Catalog Number: 177067

Product name

Rolinsatamab reference antibody

Specificity

Human PRLR / Prolactin Receptor

Antibody description

Anti-PRLR / Prolactin Receptor Antibody (Rolinsatamab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications

ELISA

Validations



SDS-PAGE

Anti-PRLR / Prolactin Receptor Antibody (rolinsatamab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.