

# ZMAT3 antibody

Catalog Number: 117067

## Product name

ZMAT3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ZMAT3 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse testis tissue, HEK-293 cells. Observed molecular weight by Western-blot: 32kd

## Preparation

This antibody was obtained by immunization of ZMAT3 recombinant protein (Accession Number: NM\_022470). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:117067(ZMAT3 antibody) at dilution of 1:800



IP Result of anti-ZMAT3 (IP:Catalog No:117067, 3ug; Detection:Catalog No:117067 1:700) with mouse testis tissue lysate 4000ug.