

# Anti-Hsp90 (4B6) Mouse antibody



Catalog Number: 168067

## Product name

Anti-Hsp90 (4B6) Mouse antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Hsp90 (4B6) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human HSP90AA1 protein fragments expressed in E.coli

## Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB

## Dilutions

WB: 1/2000

## Validations



Western blot detection of Hsp90 in HeLa, 3T3, C6, K562 and A549 cell lysates using Hsp90 mouse mAb (1:2000 diluted). Exposure time: 4min. Predicted band size: 90KDa. Observed band size: 90KDa.