

OCIAD1 antibody



Catalog Number: 113467

Product name

OCIAD1 antibody

Recommended Dilution:

WB: 1:500-1:5000

Specificity

Human, Mouse, Rat; other species not tested.

IP: 1:200-1:1000


IHC: 1:20-1:200

Antibody description

OCIAD1 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HeLa cells, human placenta tissue, human testis tissue. Positive IP detected in HeLa cells. Positive IF detected in HepG2 cells. Positive IHC detected in human liver tissue, human kidney tissue, human testis tissue. Observed molecular weight by Western-blot: 29-35 kDa

IF: 1:10-1:100

Validations




Preparation

This antibody was obtained by immunization of OCIAD1 recombinant protein (Accession Number: XM_024454108). Purification method: Antigen affinity purified.

MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:113467(ociAD1 antibody) at dilution of 1:300

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Immunohistochemical of paraffin-embedded human liver using Catalog No:113467(ociAD1 antibody) at dilution of 1:100 (under 10x lens)


Applications

ELISA, IHC, IF, IP, WB


Dilutions

OCIAD1 antibody


Catalog Number: 113467



Immunohistochemical of paraffin-embedded human liver using Catalog No:113467(ociAD1 antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using ociAD1 antibody Catalog No:113467 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-ociAD1 (IP:Catalog No:113467, 4ug; Detection:Catalog No:113467 1:300) with HeLa cells lysate 3200ug.