

# Integrin beta-7 antibody



Catalog Number: 111867

## Product name

Integrin beta-7 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Integrin beta-7 Rabbit Polyclonal antibody.  
Positive WB detected in K-562 cells, human lung tissue. Positive IHC detected in human lymphoma.  
Observed molecular weight by Western-blot: 87kd

## Preparation

This antibody was obtained by immunization of Integrin beta-7 recombinant protein (Accession Number: XM\_005268851). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:111867(ITGB7 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:111867(ITGB7 antibody) at dilution of 1:50 (under 10x lens)