

# Anti-MAGEA3 Rabbit antibody



Catalog Number: 166867

## Product name

Anti-MAGEA3 Rabbit antibody

## Specificity

Human

## Antibody description

MAGEA3 Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human MAGEA3

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of MAGEA3 in HeLa cell lysates using MAGEA3 Rabbit pAb(1:1000 diluted). Predicted band size:35kDa. Observed band size:45kDa.