

Alpha adducin antibody

Catalog Number: 107867

Product name

Alpha adducin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Alpha adducin Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive WB detected in K-562 cells, mouse testis tissue. Observed molecular weight by Western-blot: 115-120kd

Preparation

This antibody was obtained by immunization of Alpha adducin recombinant protein (Accession Number: BC013393). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107867(ADD1 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:107867(ADD1 antibody) at dilution of 1:100 (under 10x lens)