

H9N2 Neuraminidase (His Tag) recombinant protein



Catalog Number: 500167

General Information

Protein Construction

A DNA sequence encoding the influenza A virus (A/Chicken/Hong Kong/G9/97(H9N2)) neuraminidase (AAD49001.1) (His36-Ile469), termed as NA, was fused with a N-terminal polyhistidine tag.

Organism

H9N2

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant NA subunit of influenza A H9N2 (A/Chicken/Hong Kong/G9/97(H9N2)) comprises

453 amino acids and has a predicted molecular mass of 50.3 kDa. It migrates as an approximately 66-76 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H9N2 (A/Chicken/Hong Kong/G9/97) Neuraminidase / NA (His Tag) SDS-PAGE