

EBF1 antibody



Catalog Number: 110167

Product name

EBF1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

EBF1 Rabbit Polyclonal antibody. Positive WB detected in 3T3-L1 cells, Raji cells. Observed molecular weight by Western-blot: 57-67 kDa

Preparation

This antibody was obtained by immunization of EBF1 recombinant protein (Accession Number: NM_024007). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



3T3-L1 cells were subjected to SDS PAGE followed by western blot with Catalog No:110167(EBF1 Antibody) at dilution of 1:300