

# GSG1L antibody

Catalog Number: 111167

## Product name

GSG1L antibody

## Specificity

Human; other species not tested.

## Antibody description

GSG1L Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells, transfected cells. Positive IHC detected in human tonsillitis tissue. Observed molecular weight by Western-blot: 36 kDa

## Preparation

This antibody was obtained by immunization of GSG1L recombinant protein (Accession Number: NM\_144675). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



WB result of Catalog No:111167 (GSG1L antibody) with transfected cells at dilution of 1:1000.



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111167(GSG1L Antibody) at dilution of 1:700




Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111167(GSG1L Antibody) at dilution of 1:50 (under 10x lens)

# GSG1L antibody

Catalog Number: 111167

---



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111167(GSG1L Antibody) at dilution of 1:50 (under 40x lens)