

# ZNF366 antibody



Catalog Number: 117167

## Product name

ZNF366 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ZNF366 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 90-100 kDa

## Preparation

This antibody was obtained by immunization of ZNF366 recombinant protein (Accession Number: NM\_152625). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:117167(ZNF366 Antibody) at dilution of 1:300