

Mouse CD48/SLAMF2 (Fc Tag) recombinant protein



Catalog Number: 500567

General Information

Protein Construction

A DNA sequence encoding the mouse CD48 (BAE96323.1) (Met1-Arg216) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Phe 23

Molecular Mass

The recombinant mouse CD48 comprises 432 amino acids and has a predicted molecular mass of 49 kDa. The apparent molecular mass of the

protein is approximately 63.7 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse CD48 / SLAMF2 / BCM1 Protein (Fc Tag)
SDS-PAGE