

# ETV1 Polyclonal Antibody



Catalog Number: 161567

## Product name

ETV1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to ETV1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 80-300 of human ETV1 (NP\_004947.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


IHC, IF

## Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations



Immunohistochemistry - ETV1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat pancreas using ETV1 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - ETV1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using ETV1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ETV1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse pancreas using ETV1 antibody at dilution


# ETV1 Polyclonal Antibody



Catalog Number: 161567

---

of 1:100 (40x lens).



Immunofluorescence - ETV1 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using ETV1 antibody .