

HIV HIV-1 Gag-p24 (His Tag) recombinant protein



Catalog Number: 504567

General Information

Protein Construction

A DNA sequence encoding the HIV-1 gag protein P24 (Pro133-Leu363) was expressed with a polyhistidine tag at the C-terminus.

Organism

HIV

Expression Host

E. coli

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Pro 133

Molecular Mass

The recombinant HIV-1 gag protein P24 consists 242 amino acids and predicts a molecular mass of

27 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C.

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



GAG-P24 protein SDS-PAGE