

KCNIP4 antibody



Catalog Number: 111967

Product name

KCNIP4 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

human, mouse, rat; other species not tested.

Applications

ELISA, WB

Antibody description

KCNIP4 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 29kd

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of KCNIP4 recombinant protein (Accession Number: NM_025221). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111967(KCNIP4 antibody) at dilution of 1:500

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality