

ZNF256 antibody



Catalog Number: 116967

Product name

ZNF256 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZNF256 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HL-60 cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of ZNF256 recombinant protein (Accession Number: NM_005773). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116967(ZNF256 Antibody) at dilution of 1:300