

# Anti-Citrate synthetase Rabbit antibody



Catalog Number: 167967

## Product name

Anti-Citrate synthetase Rabbit antibody

## Specificity

Human, Rat

## Antibody description

Citrate synthetase Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Citrate synthetase

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, IP

## Dilutions


WB: 1/1000

IHC: 1/50

ICC/IF: 1/50

IP: 1/20

## Validations



Western blot detection of Citrate synthetase in U2OS,C6 using Citrate synthetase Rabbit mAb(1:1000 diluted)