

Anti-Phospho-c-Myc (Ser62) Rabbit antibody



Catalog Number: 167267

Product name

Anti-Phospho-c-Myc (Ser62) Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Phospho-c-Myc (Ser62) Rabbit polyclonal antibody

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser62 of human c-Myc

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000

ICC/IF: 1/200


FC: 1/20

IP: 1/20

Validations



Immunocytochemistry of c-Myc (Phospho-Ser62)(green) in HeLa cells using c-Myc (Phospho-Ser62) Rabbit pAb at dilution 1/50, and DAPI(blue)



Western blot detection of c-Myc (Phospho-Ser62) in 3T3 cell lysates using c-Myc (Phospho-Ser62) Rabbit pAb(1:1000 diluted).Predicted band size:49KDa.Observed band size:57KDa.