

# ZSCAN21 antibody



Catalog Number: 117267

## Product name

ZSCAN21 antibody

## Immunogen

[Human ZSCAN21 Recombinant protein \(GST tag\)](#)

## Specificity

human; other species not tested.

## Antibody description

ZSCAN21 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 54kd

## Preparation

This antibody was obtained by immunization of ZSCAN21 recombinant protein (Accession Number: NM\_001362779). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:117267(ZSCAN21 antibody) at dilution of 1:1000