

# MKINK2 antibody



Catalog Number: 112667

## Product name

MKINK2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MKINK2 Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive WB detected in mouse spleen tissue, MCF-7 cells. Observed molecular weight by Western-blot: 47 kDa

## Preparation

This antibody was obtained by immunization of MKINK2 recombinant protein (Accession Number: BC018345). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

## Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:112667(MKINK2 Antibody) at dilution of 1:300



Immunofluorescent analysis of (-20oc Ethanol) fixed MCF-7 cells using Catalog No:112667(MKINK2 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)