

# PDE6C antibody



Catalog Number: 113667

## Product name

PDE6C antibody

## Specificity

Human; other species not tested.

## Antibody description

PDE6C Rabbit Polyclonal antibody. Positive WB detected in Y79 cells. Observed molecular weight by Western-blot: 52 kDa; 99 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_006204). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:113667(PDE6C antibody) at dilution of 1:200